

Ceramic Repair System Kerasys LC – New syringe design starting from July, 1st

Dear Customer,

We would like to inform you that from July 1st 2019 on, we are going to change the primary package of Technovit 25... colors. From this day on, newly produced Technovit paste material will be filled in a new design syringe.

No Changes

We took care that the primary functions of the syringe are identical to the current “old” model. Even though it looks bigger, the syringe outlet has the same diameter of 3 mm.

Advantage

One big advantage of the new design is the lid of the syringe. It can be opened and closed very easily by just pulling it off and pushing it back. There is no need to screw it open anymore and still it stays on the syringe top.

Next steps

Depending on the stock situation of each Technovit 25... color, you will receive the old syringe until it is consumed or if you are going to order new colors or the stock in our warehouse is empty, you will get the new design syringe.

For any further questions, please do not hesitate to contact us again.



KULZER

MITSUI CHEMICALS GROUP

Kulzer GmbH
Leipziger Straße 2
63450 Hanau, Germany
Tel: 0049 (0)6181 9689-2574 or 2571
Fax: 0049 (0)6181 9689-3864
technik.wehrheim@kulzer-dental.com
www.kerasys.de